



Product description

Order number	55.475.001
Product description	Tube, working volume: 5 ml, (LxØ): 75 x 13 mm, material: PS, round base, transparent, push cap, natural, cap assembled, sterile, 100 piece(s)/bag

Product characteristics

Graduation	no
Cap	standard cap
Closure type	push cap
Base shape	round base
Volume of work	5 ml
Centrifugation max (RCF)	4000 x g

Size

Diameter	13 mm
Length of product	75 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Product Specification
Tube, 5 ml, (LxØ): 75 x 13 mm, PS

Material & colours

Product material	Polystyrene (PS)
Colour of product	transparent
Closure material	Low Density Polyethylene (LD-PE)
Colour of cap	natural

Purity & certification

Satisfies the requirement	sealed at 0.6 bar/10 min
Product category	no medical device IVD
Purity standard	sterile
Sterilisation	Electron irradiation
Batched	yes

Packaging

Minimum order qty.	1000
Type of smallest subpackaging	bag
Piece(s) / inner box	100
Piece (s) / outer case	1000
Piece(s) / pallet	56000
Depth of case	598 mm
Width of case	398 mm
Height of case	125 mm
Case volume	0.0298 cbm
Weight of product	0.0037 kg
Weight of case	4.35 kg
EAN of inner box	4038917107242
EAN case	4038917011723

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

