



### Product description

Order number	55.476.013
Product description	Tube, working volume: 5 ml, (LxØ): 75 x 12 mm, material: PS, round base, transparent, push cap, natural, sterile, 25 piece(s)/bag

### Product characteristics

Graduation	no
Cap	Ventilation cap
Closure type	push cap
Base shape	round base
Volume of work	5 ml
Centrifugation max (RCF)	4000 x g

### Size

Diameter	12 mm
Length of product	75 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





### **Material & colours**

Product material	Polystyrene (PS)
Colour of product	transparent
Closure material	Low Density Polyethylene (LD-PE)
Colour of cap	natural

### **Purity & certification**

Satisfies the requirement	sealed at 0.6 bar/10 min
Product category	no medical device   IVD
Purity standard	sterile
Sterilisation	Electron irradiation
Batched	yes

### **Packaging**

Minimum order qty.	1000
Type of smallest subpackaging	bag
Piece(s) / inner box	25
Piece (s) / outer case	1000
Piece(s) / pallet	48000
Depth of case	598 mm
Width of case	398 mm
Height of case	145 mm
Case volume	0.0345 cbm
Weight of product	0.0034 kg
Weight of case	4.08 kg
EAN of inner box	4038917107327
EAN case	4038917011785

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

